# DID YOU KNOW THAT THERE IS A TOOL THAT HELPS YOU EASILY FIND THE RIGHT STANDARD FOR YOUR CONNECTION?

"THANKS TO OUR MOBILE APP, DON'T HESITATE ANYMORE! EASILY DEFINE THE STANDARD WHICH CORRESPONDS TO YOUR INSTALLATION!"



#### **CONTACT US**



## Simplify your tasks

At Aseptconn, a supplier of components and products for sterile fluid transfer technology, you profit from our many years of experience in the Pharmaceutical sector.

We have developed a mobile APP that allows you to easily determine which standard is applied to your installation. Just use your phone and calliper and you will have the answer in seconds.




#### **Mobile App - Aseptconn AG**

## How to get the mobile app?

The app, which is completely free of charge, can be found in the APPLE STORE as well as on ANDROID under the name:


ASEPTCONN AG



### How does it work?

The mobile application is easy to use:

- 1. Start the application and choose your preferred language
- 2. Choose the standard for the product that you would like information about.
- 3. Measure the diameters you have available according to the image you have on the application
- 4. Make a note of the dimensions obtained from the application
- 5. Your product will correspond to the standard that is closest to the measured values



Also use the chemical compatibility table to define the seal material for your application

